



Record and analyze data with true freedom of movement using LabChart and wireless physiological monitoring.


Wireless monitoring gives your subjects unrestricted movement, allowing you to record and analyze their natural activity.

Our wireless devices and systems help you record a wide variety of signals, simultaneously. Connect seamlessly with LabChart software for streamlined data analysis.

Benefits of Wireless Physiological Monitoring:

- Record a broad range of physiological signals wirelessly via bluetooth
- Real-time streaming and analysis
- Subject comfort during recording
- Ideal for studies with single or multiple subjects
- Extended recording time with long battery life
- Suitable for close range or distance studies
- Options to record data offline then import time synchronized data for analysis


Applications include:




Autonomic




Cardiovascular



Exercise Physiology



Psychophysiology



Sleep

Typical studies:

- Apnea or Hyperapnea
- Autonomic fitness
- Autonomic function
- Breathing rate
- Cardiovascular disease and myocardial dysfunction
- Diabetic neuropathy
- Emotional reactivity
- Heart Rate Variability
- Heat Stress
- HRV / Arrhythmia
- Performance testing
- Pharmacological impact on heart health
- Physiological response
- Respiration and oxygenation during exercise



Trigno™ Wireless Foundation System

Wireless EMG in humans

Wireless EMG is ideal for capturing the intricacies of muscle movement and electrical muscle activity in subjects, especially when range of movement and comfort are important.

Whether you are studying motor control for patient rehabilitation or muscle performance, activity, and fatigue in elite athletes, ADInstruments offers streamlined wireless systems that let you record and measure electrical muscle activity simply and easily.

Direct Streaming

We offer a range of LabChart compatible solutions able to stream data directly into LabChart. With Wireless EMG studies, the Delsys Trigno Foundation system comes with both LabChart and a Trigno Device Enabler for direct data streaming.

System highlights

- Quick setup and easy to use
- Compatible with Windows
- Trigno™ Base Station holds up to 16 sensors, allowing easy future sensor additions
- Patent Pending Motion Artifact Suppression
- High resolution and sampling rate up to 2 kHz
- <500 μ s inter-sensor latency
- Real-time feedback of signal strength and battery status
- Real-time analysis options
- Wide range of analysis views and channel calculations

Foundation System Overview

RSB001DSY04

The Trigno™ Wireless Foundation System is the perfect base to build a flexible system for your movement studies.

The Trigno Base Station is equipped with 16 charge pockets which can accommodate Trigno sensors for charging and compiles data received from the active wireless sensors and transfers it over a USB 2.0 compliant connection to a Windows PC.

Contents include:

- 1 x Trigno™ Base Station Receiver (Digital)
- 1 x USB Cable
- 1 x Trigno™ Power Supply with Plug Adapter Kit
- 2 x Trigno™ Sensor Adhesive (4-slot, 90 pack)
- LabChart Pro Software
- Trigno™ Wireless Device Enabler Software



Trigno Base Station, shown with 16 Trigno sensors (purchased separately).





Trigno™ Sensor options

Select up to 16 Delsys Trigno™ sensors (sold separately). Compact and lightweight with a 40 m range, the sensors are designed for freedom of movement, allowing you to record signals directly into LabChart for analysis.



Trigno Avanti EMG + IMU Sensor

For wireless and flexible measurement of a high fidelity surface EMG signal, with a wide bandwidth 10-850 Hz and 11mV range. The gold standard for surface EMG+IMU measurements and mobile data collection.

DSY-SP-W06-14



Trigno Mini EMG + XYZ Sensor

The compact Trigno Mini Sensor is ideal for recording surface EMG on small and 'difficult-to-isolate' muscles. Applications include physical therapy, rehabilitation sciences, sports science, ergonomics, and motor control.


DSY-SP-W06-024



Trigno Snap-Lead EMG + XYZ Sensor

Connect to industry-standard 'snap on' electrode with clamp style connector leads for sEMG detection, allowing users to adjust their inter-electrode spacing as they deem suitable from difficult muscle sites.

DSY-SP-W06-018



Trigno EKG Sensor

Record a high quality ECG signal with the freedom of the Trigno system. Connects to industry standard disposable 'snap on' electrodes.

DSY-SP-W06-021



Trigno 4 Contact FSR Sensor

Record 4 independent channels of force data, each servicing an individual FSR (Force Sensitive Resistor) membrane. Ideal for recording foot pressure timing, grip force, or pressure distribution measurements.

DSY-SP-W06-020



Trigno Goniometer Adapter

For accurately measuring joint angles. Simply connect the goniometer to the adapter, activate, and begin streaming synchronized angle data.

DSY-SP-W06-023




Trigno Analog Adapter

Integrate analog outputs from dynamometers and third-party equipment. Capable of sampling and transmitting up to 4 analog channels.

DSY-SP-W06-025

Additional Accessories (purchased separately)



Trigno™ Sensor Adhesive

DSY-SC-F03



Disposable ECG Electrode packs

MLA1010/MLA1010B



Abrasive Gel

MLA1093/MLA1093B



Alcohol Swabs

MLA1094

Note: All sensors except the Analog Adapter also include a 9 Axis IMU



Equival Wireless Physiological Systems

Equival develops wearable tech products that give you accurate and precise data from real people in real environments.

Equival's compact and comfortable sensor belt and accessories are ideal for exercise research through to sleep studies. Use with LabChart's offline logging function for long sampling periods.

System highlights

- Easy setup and simple use
- High data quality with low data loss rates
- Options for live data streaming and offline data logging
- Noise and movement artifact-free ECG
- Up to 13 hours of battery life (extendable with ancillary pack)
- FDA 510(k) & CE cleared medical device

Starter Pack Overview

RSB-EQ002

Wirelessly record ECG, heart rate, expansion derived breathing rate, skin temperature and XYZ accelerometry data simultaneously into LabChart.

Contents include:

- 1 x Equival SEM
- 1 x Equival SEM USB Lead
- 1 x Equival Bluetooth Dongle
- 1 x Equival Device Enabler for LabChart (LabChart software sold separately)
- 1 x Equival Sensor Belt (RSB-EQ002)*
- or 6 x Equival Sensor Belt pack, sizes 2-7 (RSB-EQ001)

**Choose from 9 different belt sizes*



Additional Accessories (purchased separately)

- Dermal Temperature Patch **EQ-ACC-048**
- Core Temperature Capsule **EQ-ACC-023**
- Core Temperature Pill / Dermal Patch Activator **EQ-ACC-029**
- Galvanic Skin Response Sensor **EQ-ACC-034**
- Wired SpO₂ Adapter **EQ-ACC-042**
- External Battery Pack **EQ-ACC-BAT-2**
- M-Dock **EQ-ACC-MD-1**
- Additional sensor belts **EQ-02-B3**

Dermal Temperature Patch



Core Temperature Capsule



External Battery Pack



Galvanic Skin Response Sensor



Wired SpO₂ Adapter.
Connects the MLT321 SpO₂ Finger Clip to an Equival Belt.



All your data in one platform with LabChart

Use Equival's wearable tech products with LabChart software for a single, streamlined platform that lets you record and analyze multiple data sources, simultaneously.

By combining Equival products with LabChart, data can be live streamed directly into LabChart via bluetooth. Or, log recordings offline and import them at a later date for complete flexibility.



Signals

Signal	System required	Sample rate
ECG (2 Channels)	Core	256 Hz
Breathing Trace	Core	25.6 Hz
Accelerometer (3 Axis)	Core	25.6 Hz
Skin Temp (SEM)	Core	1/15 s
Skin Temp (Patch)	Core + Dermal Patch	1/15 s
Core Temp	Core + Temperature Pill	1/15 s
GSR	Core + GSR Add-On	2 Hz
SpO ₂	Core + Wired SpO ₂ Add-On	1 Hz

Typical applications

- Exercise Physiology
- Sport and Performance
- Psychophysiology
- Heart Rate Variability
- Electrocardiogram Analysis
- Sleep
- Autonomic Function

Videos and further information

Visit our blog for videos outlining:


- Configuring your SEM
- How to appropriately fit the belt to your subject
- Live streaming direct into LabChart
- LabChart's Analysis Modules and the Data Pad tool to analyze recorded data



Plus:

- Articles covering the basics of live streaming;
- How to use the LabChart logging import tool;
- A webinar exploring the inner workings of Data Pad.

For more information about our Equival range of products and solutions visit: adstruments.com/partners/equival



LabChart Data Acquisition and Analysis Software


LabChart creates a streamlined platform for all of your recording devices to work together. Acquire signals from multiple sources simultaneously. LabChart tracks every recorded action and never modifies your raw data, allowing you to easily analyze your recorded data and apply advanced calculations as your experiment unfolds.

Key Features

On	Channel Title	Sample Rate	Range	Input Settings	Units	Computed Input	Color	Style	Calculation
<input checked="" type="checkbox"/>	Channel 1	8.8	20 V	Input Amplifier...	V	Raw Data Input 1	Blue	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 2	8.8	20 V	Input Amplifier...	V	Raw Data Input 2	Red	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 3	8.8	20 V	Input Amplifier...	V	Raw Data Input 3	Green	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 4	8.8	20 V	Input Amplifier...	V	Raw Data Input 4	Yellow	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 5	8.8	20 V	Input Amplifier...	V	Raw Data Input 5	Purple	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 6	8.8	20 V	Input Amplifier...	V	Raw Data Input 6	Brown	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 7	8.8	20 V	Input Amplifier...	V	Raw Data Input 7	Pink	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 8	8.8	20 V	Input Amplifier...	V	Raw Data Input 8	Grey	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 9	8.8	20 V	Input Amplifier...	V	Raw Data Input 9	Light Blue	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 10	8.8	20 V	Input Amplifier...	V	Raw Data Input 10	Light Green	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 11	8.8	20 V	Input Amplifier...	V	Raw Data Input 11	Light Yellow	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 12	8.8	20 V	Input Amplifier...	V	Raw Data Input 12	Light Purple	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 13	8.8	20 V	Input Amplifier...	V	Raw Data Input 13	Light Brown	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 14	8.8	20 V	Input Amplifier...	V	Raw Data Input 14	Light Pink	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 15	8.8	20 V	Input Amplifier...	V	Raw Data Input 15	Light Grey	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 16	8.8	20 V	Input Amplifier...	V	Raw Data Input 16	Light Light Blue	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 17	8.8	20 V	Input Amplifier...	V	Raw Data Input 17	Light Light Green	Line	No Calculation
<input checked="" type="checkbox"/>	Channel 18	8.8	20 V	Input Amplifier...	V	Raw Data Input 18	Light Light Yellow	Line	No Calculation


Channel Settings

Get an overview of, and easily change your recording settings, calculations and channels.




Device and Channel Discovery

Easily manage signal inputs and LabChart channels all from one panel.



Units Conversion

Calibrate your recorded data to real world units.




Cyclic Measurements

Easy analysis for periodic waveforms. Find HR, systolic pressure, respiratory rate, integrate a waveform over a cycle.

A	B	C	D	E	F	G
Time	Channel 1 Avg Freq	Channel 1 Avg Rate	Channel 1 Maximum Value	Channel 1 Time at Maximum	Channel 1 Maximum Power	Channel 1 Maximum Power Frequency
1 0:01:51.0	234.8865	16589.072	131.25	0:01:51.1261	0.0	175.7813
2 0:00:00.0	302.7862	18187.1132	366.7813	0:00:00.03236	0.0002	312.5
3 0:00:01.0	254.867	15292.022	321.1563	0:00:01.1457	0.0002	312.5
4 0:00:02.0	281.7267	16603.7199	266.5625	0:00:02.2916	0.0002	312.5
5 0:00:03.0	277.3469	16640.8125	308.8125	0:00:03.2871	0.0002	312.5
6 0:00:04.0	272.2276	16333.6542	312.25	0:00:04.84415	0.0003	312.5
7 0:00:05.0	248.4031	14864.1838	276.1875	0:00:05.1304	0.0002	312.5
8 0:00:06.0	274.2485	16554.8097	328.8125	0:00:06.7861	0.0002	312.5
9 0:00:07.0	292.8983	17888.9002	307.0625	0:00:07.5086	0.0002	312.5
10 0:00:08.0	288.8575	17337.4511	275.375	0:00:08.0481	0.0002	175.7813
11 0:00:09.0	274.4335	16756.8251	337.6875	0:00:09.77968	0.0002	312.5
12 0:00:10.0	264.6496	17076.9788	210.875	0:00:10.59545	0.0002	312.5
13 0:00:11.0	276.3828	16581.1684	285.5	0:00:11.6148	0.0002	312.5


Data Pad

Analyze different segments of your time based data in a tabulated format.




Event Manager

Monitors incoming signals and detects events defined by you in order to perform a specified action.



RMS Function


Determine the energy content or intensity of your EMG signal in real time.






Specialized LabChart Modules*

*All Modules are included with LabChart Pro, or download and purchase separately.




ECG Analysis

Analyze the morphology of your ECG waveform. Automatically detects and reports the PQRST onset, amplitude, and interval in real time, or after recording. The screenshot above shows LabChart's split screen view, DVM (large numerical display of heart rate data), ECG Analysis and Cyclic Measurements.




HRV

Analyze beat-to-beat interval variation in your ECG recordings. Can be performed during recording or on a previously recorded file.



Peak Analysis

Automatic detection and analysis of multiple (non-overlapping) signal waveforms from a recording. Use in real-time or with pre-recorded data.



Video Capture

Record and synchronize a movie with your LabChart data file. Compare any data point against the video, or vice versa.

Extend your research into new territories


Please contact our expert support team to discuss a customized solution for the following applications.

Animal

Autonomic
Behavior, Sleep
and Neuroscience
Cardiovascular
Telemetry
Tissue and Circulation

In Vitro

Electrophysiology
Isolated Organ



Human

Autonomic
Cardiovascular
Exercise and Sport
Neurophysiology
Psychophysiology
Respiratory
Sleep
Speech Pathology
Tissue and Circulation

Maximize your potential

Join us for in-person training and workshops that help you learn best practice techniques and methods.

Visit adi.to/training to sign up for one of our upcoming live webinars, or access our archive of product demonstrations and application webinars on our website.

SOFTWARE TRAINING

We provide three levels of training:

- 1 The basics of data acquisition
- 2 Improving signal processing and data analysis
- 3 Automation and advanced analysis



PERSONALIZED TRAINING



APPLICATION WORKSHOPS



LIVE PRODUCT DEMONSTRATIONS



Support anywhere, anytime

Our global support system means that our expert team is always ready to help. For more information visit adi.to/support

Visit our website or contact your local ADInstruments representative for more information

ADInstruments Worldwide

Australia | Brazil | Europe | India | Japan | China | Middle East | New Zealand | North America | Pakistan | South America | South East Asia | United Kingdom

adstruments.com



ADINSTRUMENTS